


Process:
Subprocess:
Product Family:

SUSPENSIONS
OINTMENTS
BIOTRONA®



PHARMACEUTICAL

Ointment with high skin penetration efficiency



THE APPLICATION SCOPE / PURPOSE

Ointments can contain active substances that penetrate the skin when applied. The penetration rate increases with the reduction in size of this active compounds. The reduction in size often results in an increase in energy input due to high shear, which is not desired due to heat sensitivity of the compounds. Kinematica offers a solution for high turbulences for an efficient particle size reduction without high shear stresses.

Get full access to the Application Note by clicking this link:

GET FULL VERSION

Kinematica AG

Werkstrasse 7 c-d / CH-6102 Malters / Switzerland
Phone +41 41 259 65 65 / Fax +41 41 259 65 75
info@kinematica.ch / kinematica.ch

 Swiss made dispersing solutions for you. Since 1962.